

# Updating Faircom Server Licenses

## Configuration Manual





*Software by:*

PRIORITY DISPATCH CORPORATION MEDICAL PRIORITY  
CONSULTANTS, INC.


PROQA™ FOR FIRE, MEDICAL, AND POLICE

COPYRIGHT PRIORITY DISPATCH, 2007

# 1. Stop the Faircom Server service




2. Then browse to the \Program Files\FairCom\W7\W32\Server directory of the local hard drive where Faircom Server is installed and launch the Fcactvat.exe.



3. Enter in the new serial number (last six digits) and the new activation key and then click on the activate button to activate with the upgraded licenses.



4. Restart the Faircom Server service.



For additional assistance, contact Priority Dispatch Corp. Technical Support at 866-777-3911.

